

XO-101

Standing Lightweight Wheelchairs



- USA Patent / CE Approved
- 6061 T-6 Aircraft-grade Aluminum Frame
- Silver Powder Coating
- Padded Flip-back Armrests
- Flip-Up Footplates
- Clear See Through Contoured Side Panels
- Lloyd High Powered Controller Box
- Nickel Metal Hydr Battery (Included at no extra charge!)
- Easy Push Button for Stand-Up Position
- Quick Release Detachable Rear Wheels
- Thicker Corrosion/Wear Resistant Aluminum Alloy Hand Rims
- Breathable Nylon Upholstery Seat And Back w/ Pouch/Zipper Fanny Pack
- 24" x 1" High Profile Polyurethane Tires
- 8" Front Flat Free Casters
- 2-3" Memory Foam Cushion Back and Seat
- Stand-up Freeze Prevention Special Control Design
- Total Up/Down Cycles: ~250(depending on weight of user)
- Total Weight **Only: 56 lbs w/Battery**

Available In Full Power Stand-up & Propulsion Model (XO-202)!

Specifications:

XO-101N → Seat Width (16"), Seat Depth (18/19/20"), Seat Height (20"), Back Height (16-19"), Overall Height (35 – 55 ¾"), Overall Width (24 ½"), Overall Length (39 ½"), Armrest Height (9"), Weight Capacity (250 lbs)

XO-101 → Seat Width (18"), Seat Depth (18/19/20"), Seat Height (20"), Back Height (16-19"), Overall Height (35 – 55 ¾"), Overall Width (26 ½"), Overall Length (39 ½"), Armrest Height (9"), Weight Capacity (250 lbs)

*Karman Healthcare Inc. Recommends Patient Seek Physicians Counsel Before Purchasing



Karman Healthcare Inc.
19255 San Jose Ave.
City Of Industry, CA 91748-1418
Tel: 626-581-2235 Fax: 626-581-2335
www.karmanhealthcare.com Email: info001@karmanhealthcare.com